

cNIBP Protocol (v1.3)

1、Bluetooth Service Information (UUIDs):

Comm Service:	49535343-FE7D-4AE5-8FA9-9FAFD205E455
Send Characteristic:	49535343-1E4D-4BD9-BA61-23C647249616
Receive Characteristic:	49535343-8841-43F4-A8D4-ECBE34729BB3
Rename Characteristic:	00005343-0000-1000-8000-00805F9B34FB
MAC Address Characteristic:	00005344-0000-1000-8000-00805F9B34FB

PS: Host application should use the notification of the 'Send Characteristic' for package fetching.

2、Device Packet (Device to Host):

Packet length: 20 bytes

Packet Frequency: Adjustable (Default is 100Hz)

Packet Format:

Byte0: 0xff	Head mark1
Byte1: 0xaa	Head mark2
Byte2: PktIndex	Packet index
Byte3: Status	Measurement status info
Byte4: SpO2Sat	Blood oxygen saturation
Byte5: PulseRate	Pulse rate (i.e. heart rate)
Byte6: PerfuIndex	Blood perfusion index
Byte7: SBP	Systolic blood pressure
Byte8: DBP	Diastolic blood pressure
Byte9: VI	Vasoactivity Index
Byte10: SI	Signal strength Index
Byte11: PlethWave	Plethysmography wave
Byte12: Battery	Battery power in percent
Byte13: Age	Patient age
Byte14: Height	Patient height
Byte15: Weight	Patient weight
Byte16: SBP_Ref	Systolic blood pressure reference value
Byte17: DBP_Ref	Diastolic blood pressure reference value
Byte18: PacketFreq	Current packets sending frequency
Byte19: Checksum	Checksum of the whole packet

PktIndex: Range (0~255), the index of current packet, auto increment by 1

Status: 0x01 SpO2 Sensor Off

0x02 No finger

0x04 No pulse signal

SpO2Sat: Range (35~100, invalid value=127), unit: %

PulseRate:	Range (25~250, invalid value=255), unit: bpm
PerfuIndex:	Range (1~200, invalid value=0), unit: ‰
SBP:	Range (40~230, invalid value=0), unit: mmHg
DBP:	Range (40~230, invalid value=0), unit: mmHg
VI:	Range (1~100, invalid value=0), unit: %
SI:	Range (1~200, invalid value=0), unit: ‰
PlethWave:	Range (1~100, invalid value=0)
Battery:	Range (0~100), unit: %
Age:	Range (20~70), unit: years old
Height:	Range (140~190), unit: cm
Weight:	Range (40~100), unit: kg
SBP_Ref:	Range (40~230, invalid value=0), unit: mmHg
DBP_Ref:	Range (40~230, invalid value=0), unit: mmHg
PacketFreq:	Range (1, 50, 100, 200) Represent for current packets sending frequency, valid values: 1, 50, 100, 200, means 1Hz, 50Hz, 100Hz, 200Hz
Checksum	Range (0~255) Formula: Checksum = (Byte0+Byte1+...+Byte18) % 256, % is mod operator

3、Host Command (Host to Device):

Command	Byte1	description
0xff	N/A	Get software version
0xfe	N/A	Get hardware version
0xfd	Age	Set patient age value, Range (20~70)
0xfc	Height	Set patient height value, Range (140~190)
0xfb	Weight	Set patient weight value, Range (40~100)
0xfa	SBP_Ref	Set SBP reference value, Range (40~230)
0xf9	DBP_Ref	Set DBP reference value, Range (40~230)
0xf8	PacketFreq	Set packets sending frequency, Range (1, 50, 100, 200)
0xf7	on/off	Switch reference On/Off, 0x00---Off, 0x01---On

For example:

- Sending “0xff” 1 bytes, means to get software version.
- Sending “0xfe” 1 bytes, means to get hardware version.
- Sending “0xfd 0x28” 2 bytes, means set patient age value to 40 years old.
- Sending “0xfc 0xaa” 2 bytes, means set patient height value to 170 cm.
- Sending “0xfb 0x46” 2 bytes, means set patient weight value to 70 kg.
- Sending “0xfa 0x78” 2 bytes, means set SBP reference value to 120 mmHg.
- Sending “0xf9 0x50” 2 bytes, means set DBP reference value to 80 mmHg.
- Sending “0xf8 0xc8” 2 bytes, means set packets sending frequency to 200 Hz.
- Sending “0xf7 0x00” 2 bytes, means switch reference to Off.

- Software Version Package Example (20 bytes):

```
0xff 0xaa 0x53 0x56 0x31 0x2e 0x30 0x34 0x2e 0x30 0x30 0x2e 0x33 0x36 0x00 0x00 0x00
0x00 0x00 0x3a
0xff 0xaa ----- Packet head mark
0x53 ----- ASCII 'S', mark of software version
0x56 0x31 0x2e 0x30 0x34 0x2e 0x30 0x30 0x2e 0x33 0x36 ----- ASCII string
'V1.04.00.36', current software version content
0x00 0x00 0x00 0x00 0x00 ----- Padding bytes of packet
0x3a ----- Checksum byte
```

- Hardware Version Package Example:

```
0xff 0xaa 0x48 0x56 0x32 0x2e 0x30 0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00
0x00 0x00 0xd7
0xff 0xaa ----- Packet head mark
0x48 ----- ASCII 'H', mark of hardware version
0x56 0x32 0x2e 0x30 ----- ASCII string 'V2.0', current hardware version content
0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00 0x00 ----- Padding bytes of
packet
0xd7 ----- Checksum byte
```